




TAO

evolve your second nature

dentalart



sculpted, austere, intriguing
a form that appears to vanish in space





mystery surrounds a series of guarded secrets that anticipate your future needs




colour: WA7021 matt
worktop: Corian® "Glacier white"


you'll find yourself asking questions
that you've never asked before




ULTRA-SLIM FRONT & SIDES
Carved from a single 4mm thick metal sheet, drawers, door and side panels are 74% thinner than the standard 15mm. Space has also been reduced between drawers, giving a minimalist look that means more space for you and a lower total volume.



FULL EXTENSION WITH MICRO-ADJUSTMENTS
Every drawer extends completely to permit full access to content. Multiple micro-adjustments can be made in 4 directions to achieve perfect alignment.



PENTA SEAL: A STOP TO CONTAMINATION
We have used the same gaskets that we use around the frame, around the whole perimeter of the closing edge of each drawer. TAO's exclusive gaskets, similar to those found on refrigerators, are designed to collapse and adhere completely to all surfaces, ensuring an effective airtight barrier.



MAGNETIC TAGS
Say goodbye to sticky labels. TAO repositionable Tags are available in 6 colours and can be used to code drawer content or just for looks. it's up to you.

it may look familiar, but everything has changed
hidden strength

D-Box

a multi-functional tray or **SmartStorage** device?

D-Boxes can have many functions:
a drawer insert, a removable tray, or
a closed container in which to store
or transport material

D-BOX CUPS

Cups give additional organization for D-Box spaces. Available in 3 sizes and 2 heights.

D-BOX LABEL

D-Boxes can be recognized by an external colour Label that identifies its content.

D-BOX SEPARATORS



D-BOX COVER

Protects the contents during transport. D-Boxes can remain open inside a drawer as they are sealed by Penta Seal airtight gaskets.

D-BOX SEPARATORS

The 2 modular elements can enable multiple compartments – large or small – to be arranged to suit the demands of each user.

D-BOX HANDLES

More than a drawer insert, D-Box has handles that turn it into a very practical tray, easy to move so all materials are always close at hand. The handle shape offers a firm and safe hold, and, because every drawer extends completely out, extraction is simple and secure.

D-BOX LABEL

You can add an external colour Label that identifies content and is always visible even in storage cabinets or a SmartWall



D-Box

"zero tilt" extraction

simplifying actions, everyday – that's innovation

D-Box lifts vertically

THE VERTICAL GESTURE


TAO's fully opening drawers eliminate hidden zones, exposing the whole content, making it easily reachable. With D-Box, extraction is a fluid and natural as well as a "vertical gesture" that avoids obstacles or tilting.

starter kits

TAOD-BOX type 1

TAOD-BOX type 2





TAOCUPS type 1



CUP SIZES

For further spatial division, Cups are available in 3 sizes and 2 heights: type 1 for low drawers, type 2 for all other drawers.

TAOCUPS type 2



a place for everything, and everything in the
place it always should have been

Tags

say goodbye to sticky labels and say hello to magnetic Tags - the **most practical** way to identify and recognize drawer content

A identify content

TAO's magnetic Tags can serve the function of coding space for specific content

conservative hygiene endodontics



B identify users

Beyond identifying content, you can reserve dedicated spaces for different users within the same treatment room. Colours can even be associated with those used in appointment software


More than TAO's trademark look, TAO's magnetic TAGS also serve the specific function of coding space

TAO MAGNETIC TAGS These re-positionable TAGs are available in 6 colours

Labels



expand the functionality of D-Boxes by adding coloured Labels and matching them to the same colour Tag



You'll find the learning curve to colour coding rapid and intuitive, a natural way to group and store products.



Match Tags to D-Box Labels to facilitate organization of products in the treatment room

D-BOX LABELS available in the same 6 colours



the newly found beauty that you feel has
forever been part of your world





colour: red RAL 3001, grey RAL9006 (plinth)
worktop: Fenix NTM® "Bianco KDS" (white)



SmartStorage™: local

increase efficiency by creating a “**SmartSpace**” through which you can organize the passage of D-Boxes from **any** other drawers or wall cabinets

D-Box makes local storage faster, efficient, practical and protected. Every drawer, even the most distant, is now close at hand.



SMART SPACE
A transition area for placing content moved from other drawers. A central position, at the right height and distance from which you have access to what you need.



3
TO THE PLACE YOU USE THEM
Everything at hand, ready to use or protected in a closed drawer.

2
THAT MOVES FAST
Move D-Box with its entire content to the SmartSpace in a very fast and completely natural way.

1
WHOLE DRAWER CONTENT
Wherever it is, all material contained in this, or any other lower drawer, however difficult to reach, is available in an instant.

4
THEN BACK
When no longer needed, D-Box returns to its original drawer, identified by Tag and Label colour.


THE SMART REVOLUTION
The invention of D-Box has made Smart Storage™ possible. **Before** D-Box, a high-capacity container that doubled as a tray simply did not exist

SmartStorage™ : centralized


SmartWall is the perfect pairing of D-Box and centralized storage, creating a range of custom cabinets for holding a large number of D-Boxes



set up a single storeroom that supplies all treatment rooms



SmartWalls are available either open fronted or with glass door panels – transparent or smoked




Being TAO

1. metal | ultra-technical

2. glass + LED lit | a dramatic floating sensation



your new way of being. no barriers
weightless yet dynamic



3. suspended | slim profile



4. plinth | the plinth re-invented

hygiene

worktop and sink in
Corian® "Glacier White"
(X5P+LE)

worktop in
Corian® "Clay"
sink in Corian®
"Glacier White"
(X5P+LEC)

worktop in LAP
with round stainless
steel sink
(X5S+LS)

worktop in Fenix NTM®
"Bianco KDS",
with square
stainless steel sink
(XF032+LQS)

worktop in Fenix NTM®
"Grigio Londra"
with square stainless
steel sink
(XF718+LQS)

modular dispenser for disposables,
w100/150cm (WISH)

dispenser for
disposable and infra-
red soap distributor
(LUNA E L'ALTRO)

dispenser for disposables,
touch control interface (ZERO)

Living Tao

choosing an everyday
way of being. harmonies,
an irresistible feeling,
a duet of applause

selective affinities

extractable
shelf in Corian®
"Glacier White"
(X9T)

suction system and
micromotor for
touch-up work
(X13T)

keyboard drawer
(X10TT)

folding front with
extractable /
height-adjustable
aluminium arm for
trays (X330T)

D-Box module with
door. Standard or
reduced height, up
to 8 D-BOX T1/4
D-BOX T2 (TA1DB /
TA1PDB)

worktops and sinks

worktop in Fenix NTM®
"Nero Ingo", square
stainless steel sink
(XF720+LQS)

worktop and sink in
painted glass
(RV3C+RV6C)

worktop in laminate,
porcelain sink
(X5N+LP1)

handles with "matt
chrome" finish
(standard)

push&pull opening
(X108T)

accessories

drawers

mobile cabinet with
Corian® worktop
(TX6 + storage
element of choice)

Movi
an invisible technological core

search the smooth surface for imperfections
or misalignments

your experience tells you exist, but not today



rising **screen**

Your screen rises electrically only when
you need it – more table space when
don't. (fits all monitor sizes)

electric **drawers**

The wonders of a push&pull
mechanism. A simple touch opens
and closes delicately in silence

colour: grey RAL 9006
worktop: Corian® "Glacier white"

Avant system

free the space around you

flexwall™

Managing the risk of contamination of consumables

The drawer, in its lowered position, is accessible while you are seated and keeps material protected from contamination during procedures

see, check, move,

a never-ending dance of synchronized hands safely moving in a stable and harmonious universe

flexdesk™

Everything at hand

Right- or left-handed, two- or four-handed. The mobile tray follows your movements and the instruments are always there, at hand, right where you look for them.

flexclip™

Suction system no longer on chair

The ability to choose the optimal suction system position for every process allows you to maintain a better posture, facilitating your movements.

colour: white WA010 matt, grey WA7018 matt
worktop: Corian® "Glacier White"

storage cabinets

TA1

50

TA1A

25

TA1B

35

TA1D

40

TA2

50

TA2A

50

TA2B

50

TA3

50

TA3A

50

TA3B

50

TA4

50

TA4A

50

TA4B

50

TA4C

50

TA5

50

TA5A

50

TA5B

50

TA6

50

TA6A

50

TA7

50

TA14

50

TA16

50

measurements in cm

parts of TAO

sculpt it to form your inner space

hanging modules

TA1P

50

TA1AP

25

TA1BP

35

TA1CP

50

TA1DP

40

TA2P

50

TA2AP

50

TA3P

50

TA3AP

50

TA3BP

50

TA4P

50

TA5P

50

TA6P

60

TA14P

50

TA16P

50

hanging sink modules

TDA5P

70-76

TD5P

50

TD4P

50

PC cabinets

TPC3

50

TPC4

35

TB1

70-76

TB1P

70-76

sink modules

TDA5

70-76

TDA6

70-76

TD2

50

TD3

50

TD4

50

TD5

50

mobile cabinets

TX6

35

standard

hanging

avant_system

standard frame (plinth, glass / metal foot)
H69 D45cm
up to 8 drawers type1

slim frame (suspended installation)
H53 D45cm
up to 6 drawers type1

Avant frame (plinth)
H59 D45cm
up to 7 drawers type1

designing your new work place

exclusive consulting services

customization of the space

viewing the result

methods

analysis

we listen to instructions and requests, carefully assessing feasibility and relative options. Each project is unique, independently from its size, and is overseen and completed by a designer who becomes the customer's main contact person.

design and planning

items from the catalogue, made to measure or to custom specifications. Development of individual units or complete rooms, with indications for floor and wall materials and colours.


3D rendering

each project comes with a photorealistic rendering of the final result

summary: **work methods**

tailored to you, your tools, your work—
an intelligent and responsive space

the logical order of things



1

full optional
is standard


More space, safety and technology as standard. Work the way you've always worked, but better.



2

tag coding
no more sticky labels

Start using Tags to identify users and content.



3

local SmartSpace
proximity storage

Match Tags and D-Box Labels and move contents in SmartSpace drawers.



4

centralized SmartStorage™
full studio system

Store all labelled and covered up D-Boxes in a single storeroom, then distribute them to treatment rooms.

	TAO tech fundamental		organized	Advanced
THINNER STRUCTURE	●	●	●	●
DRAWER GUIDES WITH MICRO-ADJUSTMENTS	●	●	●	●
FULL EXTENSION DRAWERS	●	●	●	●
TRULY AIRTIGHT PENTASEAL	●	●	●	●
D-BOX MULTI-FUNCTIONAL INSERTS	●	●	●	●
TAGS AS CONTENT INDICATOR		●	●	●
TAGS AND LABELS			●	●
SMART STORAGE LOCAL			●	●
CENTRALIZED SMART STORAGE				●

Metal:


the irreplaceable hygiene companion

metal is an informed choice: an environmentally friendly material that provides the highest levels of hygiene

Metal is notoriously a very challenging material, but Dental Art have built themselves a reputation –in every fold and weld – for metal “designer” furniture with a distinctly “warm” feel.

WOOD exposes medical environments to the risk of cross-contamination during cleaning and maintenance. It is porous, so absorbs liquids, therefore encouraging the growth of bacterial colonies and mould on its surface and in cavities.

METAL is a non-porous material, forming a natural barrier that prevents environmental contamination, so dirt and liquid residues cannot penetrate the furniture and can be eliminated by normal cleaning operations.



specifications and customization

dimensions, materials, options, and colours

MATERIALS AND COLOURS

TAO = electro-galvanized metal sheet
D-Box cover = shockproof plexiglas
D-Box separators = shockproof plexiglas
D-Box / Cups = PST
D-Box Label = plexiglas
magnetic tags = ABS

	standard	beyond standard
	<ul style="list-style-type: none"> floor-standing: plinth floor-standing: glass / metal foot hanging: slim modules 	customized dimensions
frame	<ul style="list-style-type: none"> colour: white WAO05 	any colour glossy, matt, metallic
worktop	<ul style="list-style-type: none"> LAP, Corian®, Fenix NTM®, glass, laminate 	
drawers	<ul style="list-style-type: none"> colour: white WAO05 dimensions: from Type 1 to Type 4 soft-closing device matt chrome handles or push-pull opening system 	any colour glossy, matt, metallic electric opening system
magnetic Ttag	<ul style="list-style-type: none"> colour: grey 	orange, red, yellow, blue, green, purple
side panels	<ul style="list-style-type: none"> colour: white WAO05 	any colour glossy, matt, metallic
lower plinth	<ul style="list-style-type: none"> colour: white WAO05 	any colour glossy, matt, metallic

You can have any colour as long as it's **any colour**



Colour hues of the items in this catalogue's pictures may be slightly different from the indicated RAL/PANTONE reference, which should be considered as approximate

We believe that our furniture should blend into the environment of your studio. We have over 37,000 tones of formulated colours in stock, but feel free to add one. We adapt, you shouldn't



video TAO INTRO_Taking you somewhere new



video TAO TECH_Welcome to your workplace



AZIENDA CON
SISTEMA DI GESTIONE QUALITA'
CERTIFICATO DA DNV
ISO 9001

dentalart

made with love in Italy

dental-art.it

